

MTS Personal Care Series

# MTS ANTI-AGING COMPLEX

## 1 Product Type & Info

<b>Product Classification:</b>	MTS Facial Solution Package
<b>Product Type:</b>	MTS 0.2mm & Anti-aging products
<b>Manufactured by:</b>	Clinical Resolution Lab., Inc.
<b>Available Size:</b>	1 Box



## 2 Description

Generally, the possibility of active ingredients to penetrate the epidermal barrier (stratum corneum) through intercellular routes is approximately 0.3% per use, and due to the characteristics of cells such as keratinization, even this amount of active ingredient penetration cannot be expected all the time.

A new concept is to increase the effective penetration rate of the epidermis by forging a new transcellular route to construct a new Drug Delivery system idea.

This MTS Anti-Aging Complex provides the original and the most significant method for skin rejuvenation solution for your ultimate beauty.

## 3 Contents

1	MTS VitaPep Serum	(1 fl.oz/30ml)	Poly-Vitamin and Peptide Treatment
1	MTS Roller	(CR2: 0.2mm)	New Era of Skin Rejuvenation with 200 microneedles with GE Lexan® polycarbonate resin
5	MTS Peptide Mask		Ultimate Hydration Booster during night time
1	MTS Peptide Moisturizer	(1 oz/30g)	Regenerating, Hydrating, Smoothing Complex
1	MTS Sun Protective Cream	(1 oz/30g)	DNA Repair with SFP30
1	MTS Roller Cleaning Tray		
1	MTS Roller Cleaning Brush		
1	MTS Sterilizing Solution	(6 fl.oz/180ml)	Isopropyl Alcohol 70%, Water

## 4 Key Points

Originally Designed for Skin-Doctors.

MTS Anti-Aging Complex improves:

- Skin Texture, Tone & Pores
- Face Uplifting
- Lite Peeling
- Skin Thickness
- Whitening, Pigment Stabilizing
- Wrinkle Reduction

## 5 Direction for use

1. Step 1 - Boosting:	MTS VitaPep Serum
2. Step 2 - Correction:	MTS-Roller
3. Step 3 - Hydration:	MTS Peptide Mask / MTS Peptide Moisturizer
4. Step 4 - Protect:	MTS Sun Protective Cream
5. Step 5 - Cleaning:	MTS Sterilizing Solution